

Doc Class:	15003 - Technical Information Sheet (TIS)
Doc No:	1019473
Revision:	3
Revision Date:	24/05/2021
Created By:	4093
Document reviewed	

10000316 - Cold Water Paste Glaze Concentrate - TIS

Cold Water Paste Glaze Concentrate is a concentrated paste which, when re-constituted with hot water (Can also be made up using cold water), gives a rich, long-lasting, non-stick golden gloss when applied to all bun goods, Danish Pastries, savoury pastry lines etc.

RECIPE	kg	lb	oz	MAKE UP INSTRUCTIONS
Cold Water Paste Glaze Concentrate	1.0	2	4	Dissolve Glaze Concentrate in hot water and hand whisk until smooth. <ul style="list-style-type: none"> • Allow to cool before application. • Whisk again before use and apply evenly to defrosted pastry squares before baking, using a glaze brush. • The made-up liquid glaze should be applied by either brushing or spraying. • Ferments: Prior to, or immediately after baking. • Pastry lines: Prior to and immediately after baking. • When using a spinning disc glazing machine Cold Water Paste Glaze Concentrate can be used at a dilution of up to 7 parts water to 1 part concentrate depending on type of glaze required.
Hot/Cold Water (Hot: 60-70°C)	5.0	11	4	
Total	6	13	8	

For optimal results, using hot water gives a richer colour and shine.

Please note:

Every care has been taken to ensure that the information detailed on this sheet is correct at the time of issue. However, your recipe, other ingredients and processing are outside our control. Therefore, we advise that you carry out prudent trials to check the product makes goods of a satisfactory quality and complies with current legislation.